

REPORT DOCUMENTATION PAGE

Form Approved
OMB No. 0704-0188

Public reporting burden for this collection of information is estimated to average 1 hour per response, including the time for reviewing instructions, searching existing data sources, gathering and maintaining the data needed, and completing and reviewing this collection of information. Send comments regarding this burden estimate or any other aspect of this collection of information, including suggestions for reducing this burden to Department of Defense, Washington Headquarters Services, Directorate for Information Operations and Reports (0704-0188), 1215 Jefferson Davis Highway, Suite 1204, Arlington, VA 22202-4302. Respondents should be aware that notwithstanding any other provision of law, no person shall be subject to any penalty for failing to comply with a collection of information if it does not display a currently valid OMB control number. PLEASE DO NOT RETURN YOUR FORM TO THE ABOVE ADDRESS.

1. REPORT DATE (DD-MM-YYYY)		2. REPORT TYPE Technical Papers		3. DATES COVERED (From - To)	
4. TITLE AND SUBTITLE				5a. CONTRACT NUMBER	
				5b. GRANT NUMBER	
				5c. PROGRAM ELEMENT NUMBER	
6. AUTHOR(S)				5d. PROJECT NUMBER 2303	
				5e. TASK NUMBER M2C8	
				5f. WORK UNIT NUMBER	
7. PERFORMING ORGANIZATION NAME(S) AND ADDRESS(ES) Air Force Research Laboratory (AFMC) AFRL/PRS 5 Pollux Drive Edwards AFB CA 93524-7048				8. PERFORMING ORGANIZATION REPORT	
9. SPONSORING / MONITORING AGENCY NAME(S) AND ADDRESS(ES) Air Force Research Laboratory (AFMC) AFRL/PRS 5 Pollux Drive Edwards AFB CA 93524-7048				10. SPONSOR/MONITOR'S ACRONYM(S)	
				11. SPONSOR/MONITOR'S NUMBER(S)	
12. DISTRIBUTION / AVAILABILITY STATEMENT Approved for public release; distribution unlimited.					
13. SUPPLEMENTARY NOTES					
14. ABSTRACT					
15. SUBJECT TERMS					
16. SECURITY CLASSIFICATION OF:			17. LIMITATION OF ABSTRACT A	18. NUMBER OF PAGES	19a. NAME OF RESPONSIBLE PERSON Leilani Richardson
a. REPORT Unclassified	b. ABSTRACT Unclassified	c. THIS PAGE Unclassified			19b. TELEPHONE NUMBER (include area code) (661) 275-5015

62

separate items are enclosed

2303 MAC8

MEMORANDUM FOR PRS (Contractor/In-House Publication)

FROM: PROI (TI) (STINFO)

24 June 1999

SUBJECT: Authorization for Release of Technical Information, Control Number: AFRL-PR-ED-TP-FY99-0159
S. Tam and M.E. Fajardo, "Quantitative Matrix Isolation Spectroscopy in Heavily Doped Millimeters Thick Parahydrogen Solids"

Gordon Research Conference (International)

(Statement A)

GORDON CONFERENCE, PHYSICS AND CHEMISTRY OF MATRIX ISOLATED SPECIES
PLYMOUTH, NH, 11-15 JULY 1999

Quantitative Matrix Isolation Spectroscopy in Heavily Doped Millimeters Thick Parahydrogen Solids

Simon Tam and Mario E. Fajardo

U.S. Air Force Research Laboratory, Propulsion Directorate
(AFRL/PRSP, 10 E. Saturn Blvd. Edwards AFB CA 93524-7680

simon_tam@ple.af.mil

mario_fajardo@ple.af.mil

DISTRIBUTION STATEMENT A:
Approved for Public Release -
Distribution Unlimited

20021122 004

TASK OBJECTIVE

Develop a technique for quantifying dopant species identities and concentrations in optically dense samples using the dopant-induced infrared (IR) absorptions.

BACKGROUND

Have demonstrated we can produce gram scale samples of solid pH_2 doped with HEDM species with concentrations of 0.01 to 0.1%.

Dopants were produced using laser ablation which is not a suitable method for producing high concentrations.

Three teams in the Cryosolids Working Group tasked with: 1) developing new dopant sources; 2) developing a diagnostic for characterizing the new sources; and **3) developing diagnostic tools for detecting the products of these new sources in pH_2 .**

APPROACH

Direct absorption measurements of thick, heavily concentrated samples of HEDM doped pH_2 solids will not work as a diagnostic for these new sources.

Alternative is to use the dopant-induced IR absorptions as a diagnostic.

High Energy Density Matter (HEDM) Cryosolid Propellants

HEDM Cryosolid Program Objectives

Trap 5% molar concentration of energetic additives in solid hydrogen.
Demonstrate size-scaleable sample production method.

Payoffs

Increased Specific Impulse

$$I_{sp} \propto \sqrt{\Delta H_{sp}}$$

$$\text{LOX/LH}_2: I_{sp} = 390 \text{ s}$$

$$5\% \text{ B/H}_2 + \text{LOX}: I_{sp} = 500 \text{ s (+30\%)*}$$

*calculated for $P_{\text{chamber}} = 1000 \text{ psia}$, $P_{\text{exhaust}} = 14.7 \text{ psia}$


Greater Propellant Density

$$\text{liquid H}_2: = 0.070 \text{ g/cm}^3$$


$$\text{solid H}_2: = 0.087 \text{ g/cm}^3 (+25\%)$$

$$50/50 \text{ liquid He/solid H}_2: = 0.105 \text{ g/cm}^3 (+50\%)$$

Experimental Diagram




Ortho- and Para-Hydrogen



I.F. Silvera
Rev. Mod. Phys.,
52, 393 (1980).

B Ablation in p-H₂ Matrix Transmission Spectrum at 2 K



NOTE:
Transmission = 0.8
at 240 nm!


Demonstration of Control of o-H₂ Fraction

At 0% 0-H2:

Observation of $S_1(0)$
and non-observation
of $Q_1(0)$ implies
h.c.p. solid


Reference:

J. van Kranendonk
& H.P. Gush,
Phys. Lett., **1**, 22
(1962).




Converter Temp	% o-H ₂
135 K	70
52 K	25
37 K	8
28 K	2
15 K	<0.01


**IR Absorption Spectrum of a 37.4 μm Thick Vapor Deposited pH_2
Solid (34.4 μm front + 3.0 μm back)**




Correlation of the Total Thickness Determined by Interferometry versus the Integrated Absorption Intensity of the $Q_1(0)+S_0(0)$ Transition of Solid pH_2



Extrapolation of the Total Thickness versus the Integrated Absorption Intensity of the $Q_1(0)+S_0(0)$ Transition of Solid pH_2



Correlation between the Integrated Intensities of the $Q_1(0)+S_0(0)$ and $S_1(0)+S_0(0)$ Transitions of Vapor Deposited pH_2 Solids



NOTES

$$\alpha_{Q+S} = 188 \text{ cm}^{-2}$$

$$\alpha_{S+S} = 0.0757 \alpha_{Q+S} = 14.2 \text{ cm}^{-2}$$

For determining thickness:

Band

$Q_1(0)+S_0(0)$

Range

0-1 mm

$S_1(0)+S_0(0)$

1-10 mm

Beer's Law

$$A(\tilde{\nu}) \equiv 2.303 \log_{10} \left(\frac{I_0}{I} \right) = \alpha c l$$

$$c = \frac{A(\tilde{\nu})}{\alpha l} \Rightarrow \frac{2.303 \int_{band} \log_{10} \left(\frac{I_0}{I} \right) d\tilde{\nu}}{l \int_{band} \alpha(\tilde{\nu}) d\tilde{\nu}}$$

Increased path lengths or highly concentrated samples can cause saturation of the absorption.


If we want to work with gram scale, heavily-doped pH_2 samples, we require a spectral feature that has a very small intrinsic absorption coefficient (α) to compensate for the higher c and l .

We can use dopant-induced infrared absorptions to determine the concentration.


BUT: Need to determine $\alpha_{\text{ind}} \equiv$ the dopant-host intrinsic absorption strength

OK
Slide


Mg/pH₂ and Mg/oD₂, T=2 K



Examples of Dopant Induced H_2 Absorptions



Determining α_{ind} from α




Where:

$$\text{Slope of the line} = \frac{\alpha_{\text{ind}}}{\alpha}$$

$\alpha \equiv$ property of the dopant in the gas phase


$\alpha_{\text{ind}} \equiv$ property of the dopant and pH_2 in solid pH_2

HCl/pH₂ at 2.4 K, 88 ppm, Resolution = 0.0075 cm⁻¹
 Annealed Sample, HCl Absorptions Region




OW
5/10/03

**HCl/pH₂ at 2.4 K, 494 ppm, Resolution = 0.0075 cm⁻¹
As Deposited Sample, Induced Absorption Region**



Correlation between HCl-Induced pH_2 IR Absorption and HCl $\text{R}(0)$ Absorption



044
5/20/8

HCl-Induced pH₂ Intrinsic IR Absorption Strength

$$\int \alpha_{\text{ind}}(\text{HCl}/\text{pH}_2) dv = 0.0236 \int \alpha(\text{HCl}) dv$$

literature: $\int \alpha(\text{HCl}) dv = 19.8 \text{ km/mol}^*$

$$\therefore \int \alpha_{\text{ind}}(\text{HCl}/\text{pH}_2) dv = \underline{\underline{0.47 \text{ km/mol}}}$$

*K. N. Rao, ed., *Molecular Spectroscopy: Modern Research Vol. III* (Academic Press, Inc., New York, 1985).

Question: What is the maximum, measurable concentration of HCl/pH₂?

Assume: a) 1 mm thick sample

b) $\int A_{\text{max}} dv = 2 \text{ cm}^{-1}$

$$c_{\text{max}} = \frac{2.303 (2 \text{ cm}^{-1})}{(0.1 \text{ cm})(4.7 \times 10^4 \frac{\text{cm}}{\text{mol}})}$$

Answer:

$$= 9.9 \times 10^{-4} \text{ mol/cm}^3$$

$$\Rightarrow 2.3\% \text{ HCl}/\text{pH}_2$$

SUMMARY

For millimeters thick, heavily-doped samples, direct absorption spectroscopy fails because of limitations on dynamic range and achievable signal-to-noise levels.

Dopant-induced pH_2 transitions are a possible solution to this problem.

- 1) appear to obey Beer's Law
- 2) are very weak IR transitions (i.e., increased dynamic range for heavily doped samples)

For HCl in pH_2 , the intrinsic absorption strength of the induced transition is approximately 2.4% of the intrinsic absorption strength of HCl in the gas phase.

Can calculate the maximum measurable concentration for a HCl-doped pH_2 solid: 2.3% for a 1 mm thick sample, achieving objective of measuring $\sim 1\%$ concentration in millimeters thick samples.

FUTURE DIRECTIONS

We are in the process of completing a survey of various dopants in solid pH_2 to determine the generality of using the induced absorptions for concentration measurements.